

PowerPort^{*}

IMPLANTABLE PORT

FEEL THE NEW
STANDARD OF CARE^{*}

BARD

ACCESS SYSTEMS

Patient's Name

Doctor/Doctor's Phone

Implant Site/Implant Date

Notes

CT (300psi) **MR** (3T)

Palpation of three points on top and the triangular shape of the port help to identify the PowerPort^{*} for power injection.

Access to a PowerPort* Implantable Port must be via a PowerLoc* Safety Infusion Set if power injecting.

FLUSHING VOLUMES

- **Port not in use:**
5 ml heparinized saline, every 4 weeks
- **After each infusion of medication or TPN:**
10 ml sterile normal saline, then 5 ml heparinized saline
- **After blood withdrawal:**
20 ml sterile normal saline, then 5 ml heparinized saline
- **After Power Injection of Contrast Media:**
10ml sterile normal saline, then 5 ml heparinized saline

PowerLoc* Access Needle Flow Rates During Power Injection

Size	19 Ga.	20 Ga	22 Ga.
Color	Cream	Yellow	Black
Max Flow Rate	5ml/s	5ml/s	2ml/s
Max Pres. Setting	300 psi		

*Bard, PowerPort, PowerLoc, "Feel the New Standard of Care" and the color purple are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

Bard Access Systems, Inc.

Salt Lake City, UT 84116 USA 801-595-0700

Clinical Information Hotline: 800-443-3385

www.bardaccess.com, www.portadvantage.com

0711752 0605R

Product Code:

Lot Number:

Or place product identification sticker from the unit label here.